For your safety and continued enjoyment of this product, always read the instruction book carefully before using.

## Instruction Manual • Model DP18

DIRECT PLUG-IN TOOTHBRUSH

by CONAIR®

# INTERPLAK®

## Using the Interplak® Direct Plug-In Toothbrush: The Fastest and Easiest Way to Brush

### 6 Simple Steps:

#### 1. Charge the lower battery unit.


**IMPORTANT NOTE:** The plug blades and face must be clean and dry before plugging into AC outlet to reduce risk of burns, electrocution, fire or injury. Allow base to charge for at least **24 hours** before first use.

**Note:** If the outlet is controlled by a wall switch, be sure the switch is always on for a continuous charge. Plug the lower battery unit into a live outlet and leave it charging for 24 hours. Use only a 120v, 60Hz outlet.

#### 2. After recharging:

Take lower battery unit out of outlet and line up prongs on base to the opening in the handle and push together until in place. The plug blades are polarized (one blade is wide) and can go into the handle only one way. Push together until there is no gap between handle and base.

**IMPORTANT NOTE:** The unit will not function with a slight gap between the handle and the base. Make sure that the base is fully pushed in.

#### 3. Attach the brush head.

Attach and remove with just a simple twist.

#### 4. Wet the bristles, apply a pea-size amount of toothpaste, and put in mouth.

Clear gel or low abrasion formulations are strongly recommended to enhance the life and performance of your brush head, but pastes may also be used. To avoid splashing, switch the power on after placing the brush head in your mouth.

Conair will (at our option) repair or replace your appliance free of charge for 24 months from the date of purchase if the appliance is defective in workmanship or materials. To obtain service under this warranty, return the defective product to the service center listed below, together with your sales slip and \$5.00 for postage and handling. California residents need only provide proof of purchase and should call 1-800-326-6247 for shipping instructions. In the absence of a purchase receipt, the warranty period will be 24 months from the date of manufacture.



ANY IMPLIED WARRANTIES, OBLIGATIONS OR LIABILITIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SHALL BE LIMITED IN DURATION TO THE 24 MONTH DURATION OF THIS WARRANTY. Some states do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you.

IN NO EVENT SHALL CONAIR BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS OR IMPLIED, WHATSOEVER. Some states do not allow the exclusion or limitation of special, incidental or consequential damages, so the above limitation may not apply to you. This warranty gives you specific legal rights, and you may also have other rights which vary from state to state. For service, please bring or send defective unit with sales slip to the nearest service center listed below.

**Limited 2 Year Warranty**

**5. Turn on and go slowly.**

Press the button to turn on the unit. Press it again to turn it off. Most dentists recommend brushing for 2 minutes. Place the bristles perpendicular to flat surface of tooth, and **slowly** guide the brush from tooth to tooth, pausing as you go so bristles can work into crevices. The Interplak® Direct Plug-In Toothbrush will automatically provide the correct cleaning action so there is no need to scrub back and forth. Be sure to run brush along your gumline. Brush all surfaces, including the outsides, insides, and chewing surfaces. Pay particular attention to back molars.

**6. IMPORTANT:**


**Rinse the brush head.**

After each use, leave the unit running, and run warm water through the cleaning ports on front and back for 15 seconds to thoroughly rinse out toothpaste residue.

**VERY IMPORTANT:** If not cleaned properly, toothpaste may harden inside the brush head and reduce its life. Towel-dry power handle and brush head to keep clean and free of build-up.

**IMPORTANT NOTE:** Do not take unit apart when rinsing.

Do not rinse plug blades or face. Plug blades and face must be clean and dry before plugging into AC outlet.

**Replace your brush head at least every 3 to 4 months.**

Over time, the gears in the brush head become worn, reducing the power and effectiveness of the unit. You may notice the unit slowing down, making a clacking noise, or jamming.

# OPERATING INSTRUCTIONS

**MEDICAL WARNING** For your safety:

- Contact your dentist prior to use if you have had oral or periodontal surgery in the previous twelve months.
- When you start any new oral hygiene program, your gums may initially bleed slightly. Contact your dentist if excessive bleeding occurs or if bleeding continues after two weeks of use.
- Be sure that only the bristles touch your teeth and gums. Damage to teeth and gums may occur if any other part of the unit touches them.
- This product is designed for oral cleaning only.
- DO NOT use the unit for any other purpose.
- As with any other new oral hygiene program, individuals with heart problems should check with their physician prior to using the INTERPLAK® Direct Plug-In Toothbrush, and should always exercise particular care with any oral hygiene device.

**WARNING**  
Changes or modifications to this unit not expressly approved by Conair could void the user's authority to operate the equipment and void limited warranty. Servicing should always be referred to an Authorized Service Center.

- Avoid using a brush with crushed or bent bristles. Damaged bristles may break off while the unit is in use. In order to prevent this, replace your brush head if bristles are worn or bent. Dental professionals recommend replacing your brush head every three to four months or sooner if the bristles become worn or bent.
- The Interplak® Direct Plug-In Toothbrush is a personal care unit and is intended for household use only.

## Common Answers and Answers

### My gums are tender and bleed a little. Is this normal?

During the first few days of using the Interplak® unit, your gums may bleed slightly. This is common when you begin a new oral hygiene program. After a few days, the bleeding will generally cease, and you will be on your way to healthier gums. If your gums continue to bleed after two weeks, see your dentist.

### How often should I use my Interplak® Direct Plug-In Toothbrush?

Use the Interplak® Direct Plug-In Toothbrush to replace your manual brushing routine, twice daily or as recommended by your dentist.

### How do I clean the unit?

Unplug the base from the outlet, and clean with a cloth or tissue moistened with warm water. Submerging the unit or using a dishwasher is not recommended and may cause damage or injury.

**IMPORTANT NOTE:** The lower battery unit must be clean and dry before plugging into AC outlet.

### Should I run down the battery periodically?

The Interplak® unit has a nickel-cadmium (Ni-Cd) battery similar to those in many cellular phones and electric razors. To maintain the full charging capacity and extend the life of the battery, it must be periodically run down during normal use until it starts to slow down, then recharged.

### Why are there different colored brush heads?

The brush heads for your model are identical except for the color indicator ring. This helps you keep track of which family member is using which brush head.

### Whom do I call for questions?

You may contact one of our customer service representatives at 1-800-334-4031 between 8 AM and 5 PM Mountain Standard Time, Monday through Friday.

# SAVE THESE INSTRUCTIONS

## WARNING

To reduce the risk of burns, electrocution, fire, or injury:

- Never operate the unit if it has a damaged plug, if it is not working properly, if it has been dropped or damaged, or if it is immersed in water. Refer to an Authorized Service Center for repair service.
- Unplug the lower battery unit by grasping the unit firmly.
- This product should not be used by, on, or near children or individuals with certain disabilities.
- Use this product only for its intended use as described in this manual. DO NOT use attachments unless recommended by the manufacturer.
- Never use while drowsy.
- DO NOT use the unit while bathing.
- DO NOT place or store the lower battery unit where it can fall or be pulled into a tub, sink, or toilet.
- DO NOT use after the lower battery unit has fallen into water.
- DO NOT submerge the lower battery unit in water or any other liquid.

**DANGER**  
To reduce the risk of electrocution:

## READ ALL INSTRUCTIONS BEFORE USING IMPORTANT SAFEGUARDS

### Troubleshooting

#### **The instrument will not turn on.**

- Take the brush head off the handle and then turn the unit on. If the shaft on the handle is rotating rapidly, the unit is still operational, but you will need to replace your brush head.
- If the shaft is not rotating, push the lower battery unit tightly together so there is no gap. Then turn the unit on.
- If the shaft on the handle is not rotating rapidly, your battery may be discharged. Check to see if the charger is plugged into a live outlet. If plug is controlled by a wall switch, make sure it is switched on during charging. The LED light should illuminate. Charge the unit for 24 hours and try again.
- If still unsuccessful, lower battery unit may need to be replaced.

#### **The bristles stop rotating during use.**

- You may be applying too much pressure against the teeth. Unlike a manual brush, the Interplak® Direct Plug-In Toothbrush cleans effectively without strong pressure.
- Try cleaning the brush head thoroughly with warm water or mouthwash.

#### **My brush head stopped working.**

- If your brush head stopped working soon after you began using it, try using a different toothpaste and rinse thoroughly after each use. Regular gel toothpastes provide excellent cleaning and extend the life of brush heads.

#### **The power handle is warm to the touch.**

- This is normal, caused by the system used to charge the instrument's batteries.

### For Battery Disposal

THIS PRODUCT HAS ENCLOSED BATTERIES IN THE LOWER BATTERY UNIT. PLEASE RECYCLE OLD BATTERIES TO HELP PROTECT THE ENVIRONMENT.

#### **Contains nickel-cadmium battery**

Nickel-Cadmium batteries should not be put in municipal waste systems. Contact your local waste collection agency for proper recycling and disposal information for your area.

Since this product has no user-serviceable parts, for lower battery unit replacement call 1-800-334-4031 between 8 AM and 5 PM Mountain Standard Time, Monday through Friday.

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Conair Corporation